


Triple Layered Disposable Face Masks

Nyeleti Health Catalogue	Box of 50, excl VAT, delivered	Box of 50, incl VAT, delivered
Pack of 50 Disposable Surgical Face Masks (Medical): <ul style="list-style-type: none"> • 3 Ply Suitable for use in Health Facilities and High Infection Environments	R 527.63	R 606.77
Pack of 50 Disposable Face Masks (Non Medical): <ul style="list-style-type: none"> • 3 Ply, • Dawn Type DN01. Suitable for use in Construction, Manufacturing, Schools, Offices, Security Services, Transport Trains, Buses and Taxis	R 127.86	R 147.04



Type: DN01

Application: nose and mouth protection.



Outer layer: non-woven isolation

Filter layer: electrostatic adsorption melt-blown filter

Inner layer: skin friendly non-woven fabric

Executive standard: GB/T 32610-2016

Product specification: 175mm×95mm (±5mm), 50pcs/box

Applicable: adults

Shelf life: 2 years



Usage

1. Spread out the mask, put the elastics on your ears.
2. Fix the nose clip in the shape of the nose.
3. Please pull the mask to the jaw.

Warning

1. Should any damage be identified, DO NOT USE IT.
 2. Failure to follow the usages of this product could reduce the effectiveness of the respirator.
 3. It is strictly forbidden to use it overdue.
- Storage conditions: stored in a cool and dry place, humidity less than 80%, keep away from fire and pollution.